

# A Fast Protocol for the Extraction of High Molecular Weight Nuclear DNA (>50 kbp) from Plants

**Note:** This method is combined method of nucleus extraction plus CTAB method. It is cp- and mt-genome depleted.

**Note:** This method is not recommended for samples which produce a viscous solution when ground and mixed with isolation buffer.

## I. PREPARATION OF SOLUTIONS

### 1. Preparation of nuclei isolation buffer (IB) (for 10 reactions)

- For 200 ml of nuclei isolation buffer, dissolve the following in ca. 100 ml water

Tris-HCl (1 M stock; ph 9.5) 3 ml (Final conc.: 15 mM)

EDTA (0.5 M stock) 4 ml (Final conc.: 10 mM)

KCl 1.94 g (Final conc.: 130 mM)

NaCl (5 M stock) 0.8 ml (Final conc.: 20 mM)

- Add 16 g of PVP-10 gradually while spinning the solution fast with a magnetic stir bar.

- Fill with water to 200 ml (enough for 5 extractions).

- Add 0.05 g of spermine and 0.07 g of spermidine (Store IB at 4 °C).

- After mixing IB with ground sample,

add 0.1 % Triton X-100 and 7.5 % B-mercarptoethanol to IB (this constitutes IBTB)

(e. g., to 40 ml of IB buffer, add 40 µl of Triton X-100 and 3 ml of B-mercarptoethanol).

**Note:** Keep in 4C refrigerator and use it within two weeks will be safe.

### 2. Preparation of Carlson lysis buffer

- Carlson lysis Buffer = 2X CTAB (cetyl trimethylammonium bromide) buffer + 1% PEG 6000

- For 100 ml of Carlson lysis buffer

Tris-HCl (1 M stock ph 9.5) 10 ml (Final conc.: 100 mM)

EDTA (0.5M stock) 4 ml (Final conc.: 20 mM)

NaCl 8.2 g (Final conc.: 1.4 M)

CTAB 2 g (Final conc.: 2 %)

PEG 1 g (Final conc.: 1 %)

**Note:** Keep in 4C refrigerator and use it within a month will be safe.

## II. GRINDING AND NUCLEI ISOLATION

1. Cool down mortar and pestle in -80 °C before the experiment starts. Grind 2g of tissue to a fine powder under liquid Nitrogen.

2. Add 2 g of ground leaves into 20 ml IB buffer in a 50 ml conical tube and mix by inverting.

**Note:** Grinding level (fine/coarse) doesn't seem to appreciably affect to HMW DNA size but insufficient grinding appears to reduce the extraction efficiency and DNA yield.

**Note:** Increase amount of ground sample in case of succulent plants

3. Immediately add 20 µl of Triton X-100 and 1.5 ml of B-mercarptoethanol and mix by inverting.

4. Keep on ice for 10 min.

**Note:** this step should be conducted inside a fume hood because the IBTB contains B-mercarptoethanol, which is toxic.

5. Filter mixture through a 100  $\mu\text{m}$  vacuume-aid cell strainer seated on a 50  $\text{ml}$  conical tube to collect nuclei suspension (Fig. 1). To aid filtration, gently scrap unfiltered plant tissue on the filter with top side of 1000  $\mu\text{l}$  pipette (blue) tip. The filtrate should look light green color.
6. Repeat the filtering step again with 40  $\mu\text{m}$  cell strainer.
7. Add 200  $\mu\text{l}$  Triton X-100.  
**Note:** This is a step to lyse cell and organelle membranes except nuclear membrane.
8. Centrifuge for 10 min. with 3000g at 4  $^{\circ}\text{C}$ . to pellet the nuclei (Fig. 2).
9. Discard the supernatant.

### **III. EXTRACT NUCLEAR DNA USING CTAB buffer**

1. Add 5  $\text{ml}$  of Carlson Lysis buffer + 12.5  $\mu\text{l}$  B-mercaptoethanol to the tube and re-suspend nuclei pellet with brief tapping.  
**Note:** Incomplete resuspension could reduce yield since many nuclei are not lysed by CTAB. Pipette pellet briefly with end-cut 1,000  $\mu\text{l}$  pipette tip to help resuspension and vortex sample.
2. Incubate at 65  $^{\circ}\text{C}$  for 15 min. (maximum 2 hours)  
**Note:** If it is not resuspended after the incubation, centrifuge it (3000g for 5 min.) and just use supernatant for faster processing.
3. Transfer suspended nuclei pellet to 15  $\text{ml}$  polypropylene tube and add 5  $\text{ml}$  (equal volume) of chloroform : isoamyl alcohol (24 : 1)
4. Invert several times to mix.
5. Centrifuge (3000 g) for 10 min. at 4  $^{\circ}\text{C}$ .
6. Transfer aqueous upper phase to a new tube.  
**Note:** Take only 9/10 of supernatant. This step is highly correlated to the quality of extracted DNA.  
**Note:** If it is viscous, slow pipetting will help not to suck up plant tissue.
7. Repeat 3~6 (optional but highly recommended).
8. Add 1/10 volume of 3 M NaOAc, mix gently, add some volume of room temperature isopropanol, and gently invert several times.  
**Note:** ex) For 4.5  $\text{ml}$  of supernatant, add 0.45  $\text{ml}$  of 3 M NaOAc and add 4.95  $\text{ml}$  of isopropanol.
9. Precipitate at -20  $^{\circ}\text{C}$  more than an hour.  
**Note:** If you can see precipitate, go to the next step for the fast processing.
10. Centrifuge (3,000 g) for 10 min. at 4  $^{\circ}\text{C}$ .
11. Discard supernatant.
12. Wash pellets with 70% cold ethanol (ca. 20 ml per a tube).
13. Centrifuge (3,000 g) for 10 min. at 4  $^{\circ}\text{C}$ .
14. Discard supernatant.  
**Note:** keep tube inverted for a min. and wipe out tube wall with Kimwipe.
15. Air dry pellet completely.  
**Note:** this step is very important for the quality of DNA.

#### **IV. RNase A and Protease K treatment**

1. Dissolve pellet with 2 ml of TE buffer

**Note:** If it is difficult to dissolve, put the tube in 50 °C water bath for a while.

**Note:** Crushing pellet with pipette tip might be helpful for fast dissolving but never vortex sample. If it is not dissolved after the incubation, centrifuge it (3000g for 5 min.) and just use supernatant for the fast processing.

2. Add 20 µl (10 ul/ml) of RNase A (10 mg/ml conc.)
3. Incubate 37 °C 5 min.
4. Add 20 µl (10 ul/ml) of proteinase K (>600 units/ml conc.)
5. Incubate 50 °C 15 min.
6. Add 2 ml (equal volume) of chloroform : isoamyl alcohol (24 : 1)
7. Invert several times to mix
8. Centrifuge (3000 g) for 10 min. at 4 °C.
9. Transfer aqueous upper phase to a new 15 ml tube .

**Note:** Take only 9/10 of supernatant. This step is highly correlated to the quality of extracted DNA.

10. Repeat 6~9 (optional). Use 2 ml Eppendorf tube in the final repeat.
11. Add 1/10 volume of 3 M NaOAc, mix gently, add some volume of room temperature isopropanol, and gently invert several times.

**Note:** ex) For 1.7 ml of supernatant, add 0.17 ml of 3 M NaOAc and add 1.87 ml of isopropanol.

12. Precipitate at -20 °C more than an hour.

**Note:** If you can see precipitate, go to the next step for the fast processing.

14. Centrifuge (3000 g) for 10 min. at 4 °C.
15. Discard supernatant.
16. Wash pellets with 70% cold ethanol (ca. 5 ml per a tube).
17. Centrifuge (3000 g) for 10 min. at 4 °C.
18. Discard supernatant.

**Note:** keep tube inverted for a min. and wipe out tube wall with Kimwipe.

#### **19. Completely drying.**

**Note:** this step is a critical step for the quality of DNA.

20. Add 50~ 500 µl of water in each tube to dissolve pellet.

**Note:** If it is difficult to dissolve, put the tube in 50 °C water bath for a while.

**Note:** Crushing pellet with pipette tip might be helpful for fast dissolving but never vortex sample. If it is not dissolved after the incubation, centrifuge it (3,000g for 5 min.) and just use supernatant for the fast processing.

21. Check quality (A260/280 and A260/230) and quantity of extracted DNA using Nanodrop and Qubit.
22. Check length distribution using Femto pulse.

## **ORDER**

### **- PVP-10**

Sigma-aldrich CAS Number: **9003-39-8** 100G \$43.2

<https://www.sigmaaldrich.com/catalog/search?term=PVP-10&interface=All&N=0&mode=match%20partialmax&lang=en&region=US&focus=product>

### **- spermine**

Sigma-aldrich **S2876-1G** \$47.5

<https://www.sigmaaldrich.com/catalog/search?term=S2876&interface=All&N=0&mode=match%20partialmax&lang=en&region=US&focus=product>

### **- spermidine**

Sigma-aldrich **S2501-1G** \$31

<https://www.sigmaaldrich.com/catalog/search?term=S2501&interface=All&N=0&mode=match%20partialmax&lang=en&region=US&focus=product>

### **- Cell strainer**

Sigma-aldrich **Z742101 yellow 100 µm** pk of 50 \$94.5

<https://www.sigmaaldrich.com/catalog/product/sigma/z742101?lang=en&region=US>

### **- Cell strainer**

Sigma-aldrich **Z742102 blue 40 µm** pk of 50 \$94.5

<https://www.sigmaaldrich.com/catalog/product/sigma/z742102?lang=en&region=US>

### **- PEG 6000**

Sigma-aldrich 81260-1KG \$68

<https://www.sigmaaldrich.com/catalog/product/aldrich/81260?lang=en&region=US>

### **- RNase A**

Sigma-aldrich **R6513-10MG** \$35.75

### **- Proteinase K**

Sigma-aldrich **P2308-5MG** \$43

<https://www.sigmaaldrich.com/catalog/search?term=Proteinase+K&interface=All&N=0&mode=match%20partialmax&lang=en&region=US&focus=product>